Alivion Spark M-30

EN

Instructions for use
📄 1

DE

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FR

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ES

Manual de uso
📄 87

IT

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Most up-to-date version of the manual
and translation in other languages
available online at

www.alivion.ch.

Contents

- 1 For your safety..... 2**
 - 1.1 General safety instructions..... 2
 - 1.2 Battery..... 3
 - 1.3 Electrical safety 3
 - 1.4 Attention signs 4
 - 1.5 Liability 4
- 2 Description..... 5**
 - 2.1 Product overview 5
 - 2.2 Intended use..... 6
 - 2.3 Required consumables 6
- 3 Use..... 7**
 - 3.1 Recharging the battery 7
 - 3.2 Switching the device on and off 8
 - 3.3 Navigation and main menu 9
 - 3.4 Measurement 10
 - 3.5 Calibration 13
 - 3.6 Sugar correction..... 14
 - 3.7 Warnings..... 15
- 4 Data..... 16**
 - 4.1 On the device 16
 - 4.2 Export..... 17
- 5 Settings..... 18**
- 6 Device faults 19**
- 7 Troubleshooting 20**
 - 7.1 Sampling errors..... 20
 - 7.2 Measurement errors 21
 - 7.3 Calibration errors..... 22
 - 7.4 General errors..... 23
- 8 Maintenance..... 25**
 - 8.1 Transport and storage..... 25
 - 8.2 Cleaning 25
 - 8.3 Disposal..... 26
- 9 Technical data 27**

1 For your safety

1.1 General safety instructions

- Before using the product, carefully read these instructions for use and those of the associated products.
- Follow the instructions for use carefully. The user must fully understand these instructions and follow them exactly. The product may only be used for its intended purpose.
- Do not discard the instructions for use. Ensure storage and proper use by the user.
- Follow local and national regulations that apply to this product.
- Only personnel trained by Alivion may check, repair and maintain the product.
- Do not use defective or incomplete products. Do not make any changes to the product.
- Make sure that the sensor is not dropped or subjected to other strong vibrations. As a result, the sensor could be damaged without any visual damage.
- Inform Alivion in the event of unexpected errors or failures (not listed in chapter 7) of the product or product parts.
- Only use original Alivion parts and accessories for maintenance work. Otherwise, the correct functioning of the product could be impaired.
- The device is not approved or certified for operation in explosive atmospheres.
- Electrical coupling with devices that are not mentioned in these operating instructions should only be carried out after consulting Alivion.

1.2 Battery

- The device contains a rechargeable lithium-ion polymer battery. Do not replace the batteries yourself. The batteries must be removed and replaced by a qualified professional. Please contact Alivion's technical support for more information.
- Stop using the device if a battery is leaking. If battery fluid comes into contact with your eyes, immediately rinse your eyes with plenty of water and consult a doctor. If any fluid leaks on your hands, wash them thoroughly with water. Carefully wipe the fluid from the exterior of the device with a cloth.

1.3 Electrical safety

- Only use the provided USB charger to charge the device.
- Do not expose the device to fire, microwaves, high temperatures or direct sunlight.
- Do not let the device come into contact with liquid and do not use it with wet or oily hands. If liquid gets inside, stop using and contact 's technical support.
- If you hear a strange noise, see smoke or smell something strange, unplug the USB charger from the socket, cancel running measurements and turn off the device. Do not turn on the device again and contact 's technical support.

1.4 Attention signs

The following warning signs are used in this document to identify and emphasize related warning texts that require special attention by the user. The meanings of the warning signs are defined as follows:

WARNING

Indication of a potentially dangerous situation. If this is not avoided, injuries or damage to the product or the environment can occur. Can also be used as a warning against improper use.

NOTICE

Additional information on the use of the product.

1.5 Liability


The liability for the proper function of the device is irrevocably transferred to the owner or operator to the extent that the device is improperly serviced or repaired by personnel not employed or authorized by Alivion or if the instrument is used in a manner not conforming to its intended use.

Alivion cannot be held responsible for damage caused by non-compliance with the recommendations given above.

The warranty and liability provisions of the Alivion's General Terms and Conditions of Sale are likewise not modified by the recommendations given above.

2 Description

2.1 Product overview

Content delivered with the device

- A. USB charger
- B. USB cable
- C. Spare sealing O-ring

2.2 Intended use

The Alivion *Spark M-30* is used to quickly measure the **methanol concentration** in aqueous liquids with **ethanol** concentrations up to **20 vol%**. This includes soft drinks, wines & beers and their mixtures with soft drinks, fruit juices, etc. The device is handheld, suitable for mobile and network-independent operation.

2.3 Required consumables

The following consumables are required for proper operation:

- Alivion Methanol Standards 1 and 2 WITHOUT ETHANOL
- Alivion Methanol Standards 1 and 2 WITH ETHANOL
- Alivion liquid sample vials and pipettes

These products may be included already in your starter kit and can be ordered separately from Alivion.

WARNING

Accurate results can only be ensured if official consumables purchased by Alivion are used.


3 Use

- The device is delivered with a factory calibration. However, Alivion strongly recommends performing a **calibration** before the **first measurement** in order to achieve the maximum measurement accuracy.
- The device must be **calibrated regularly**. Alivion recommends performing a new calibration at the beginning of **each measurement day** in order to guarantee maximum accuracy.
- **Carbonated** samples can be measured directly, however slightly lower measurement accuracy must be expected. For best results, **de-carbonate samples** before measurement.
- The device can be used for measurements within the operational temperature range of **10–35 °C** and delivers accurate results using **temperature correction** algorithms.
- While the device accounts for temperature variations, the measurement is **highly temperature sensitive**. Please carefully follow the instructions in **chapter 3.4** to achieve maximum accuracy.

3.1 Recharging the battery

Charge the device by connecting it to a wall socket using the provided **Alivion USB charger**.

The device battery level is indicated with the corresponding icon in the top right part of the screen.

NOTICE

Measurements are not possible during charging.

- The Li-ion polymer battery can be charged at any time. Interrupting the charging process will not damage the battery.

- Avoid keeping the battery at 100% all the time, as this may cause accelerated aging. Therefore, it is recommended to charge the battery if the charge level is $\frac{1}{3}$ or less.
 - When stored for a long time, the battery should neither be empty nor fully charged. Be sure to charge the built-in battery at least every six months. If the battery is not used for a long period of time, it may not be able to be charged.
 - Using or storing the Spark M-30 at temperatures above 35°C or below 0°C may accelerate the aging of the battery.
-

**NOTICE**

The battery is shipped partially charged. It is recommended to fully charge the battery with the supplied charger before first use.

3.2 Switching the device on and off

- To switch the device on and off, press and hold the **home button** for 3 seconds.
 - When turned on, please wait a few seconds for the sensor to heat up and prepare.
 - The warning flag **Tem** typically appears in the first minutes after starting the device due to internal heat-up. Refer to chapter 3.7 for more details.
-

**NOTICE**

The noise coming from the device is from the internal pump and it is normal.

3.3 Navigation and main menu

The user interface is navigated by the 4 face buttons:

Confirmation

Confirms the current highlighted selection.


Navigation

Change the current selection (left/right or up/down).

Home

Shortcut back to the main menu. Also turns the device on/off when held for 3 seconds.

3.4 Measurement

A measurement is started by selecting the option *Measure* in the main menu and following the instructions that appear on the display. A measurement always consists of:

1. Sample preparation
2. Sampling
3. Analysis
4. Result
5. Regeneration


Sample preparation

Always prepare samples freshly right before measurements.

Your sample liquid, the vials and pipettes **must be at ambient temperature** to avoid measurement errors.

Reusing sample vials or pipettes will lead to wrong results.

1. Withdraw **2 mL** of your sample liquid with a fresh pipette.
2. Release the liquid into the vial and onto the filler material.
3. Wipe the rim of the sampling vial with a clean tissue or your hand.

WARNING

Do not wet the windings or rim of the vial. This will lead to inaccurate results and can poison the device. If the device comes into contact with sample liquid, stop the measurement, clean the areas with a wet cloth and let the device dry before the next measurement.

Sampling

For sampling, follow the instructions displayed on the screen.

1. Screw the sample vial into the bottom of the device.
2. Wait for sample extraction.
A timer will show the progress and a sound is played when extraction is completed.
3. Unscrew the sample vial and confirm.

WARNING

Proceed with sampling directly after sample preparation.

WARNING

Avoid excessive touching of the sample vials to prevent heat-up.

Analysis

Analysis usually takes between 2 to 4 minutes, depending on the ambient temperature. A sound is played and the result is displayed as soon as the analysis is finished.

Results

At the end of the analysis, the results screen appears and the measurement can be saved or deleted.

- Every saved measurement is given a sequential ID number.
- Saved measurements can be re-viewed from the *Data* menu (see chapter 4.1) and exported to a PC (see chapter 4.2).

The screenshot shows a measurement screen with the following data and callouts:

- Measurement details:** Points to the 'Nr.', 'Date', and 'Time' fields.
- Methanol concentration:** Points to the large '240' value.
- Temperature:** Points to the 'Temp, °C' field showing '25.2'.
- Sugar concentration (chapter 3.6):** Points to the 'Sug, g/L' field showing '150'.
- Warnings (chapter 3.7):** Points to the 'Warnings: Tem' indicator.
- Details to corresponding calibration:** Points to the 'Calibration' section showing 'Date', 'Time', and 'Temp, °C'.

At the bottom, there are two buttons: 'DELETE' (with a trash icon) and 'SAVE' (with a document icon).

NOTICE

The methanol concentration is displayed in parts per million, ppm (mg/L).*

Regeneration

The device must regenerate between each measurement. It is done automatically and regeneration time varies between several seconds to several minutes, depending on the ambient temperature and the measured sample. Please note that regeneration takes significantly longer with higher ethanol content of the measured sample.

NOTICE

If the device is turned off during regeneration, it will be completed when starting the device the next time.

* Mass concentration (DIN 1310) at a reference temperature of 25 °C.

3.5 Calibration

A calibration is performed similarly to the sample analysis described before in chapter 3.4 using the *Alivion Methanol Standard 1* and 2 *WITHOUT/WITH ETHANOL*.

- A measurement is always evaluated with the most recent saved calibration.
- Regular calibration is essential to ensure maximum accuracy.
- If a new calibration should be performed, the warning **Ca1** is displayed on the home screen (see chapter 3.7 for details).

To perform a calibration, select *Calibrate* from the main menu followed by the calibration mode:

Calibration without ethanol

- Measurement with *Alivion Methanol Standard 1 WITHOUT ETHANOL* (orange liquid, red bottle cap)
- Measurement with *Alivion Methanol Standard 2 WITHOUT ETHANOL* (blue liquid, red bottle cap)

Gives the maximum measurement accuracy when measuring samples that do not contain ethanol.

Calibration with ethanol

- Measurement with *Alivion Methanol Standard 1 WITH ETHANOL* (orange liquid, blue bottle cap)
- Measurement with *Alivion Methanol Standard 2 WITH ETHANOL* (blue liquid, blue bottle cap)

Gives the maximum measurement accuracy when measuring samples that contain ethanol.

NOTICE

Alivion recommends performing a calibration in the beginning of each measurement day to achieve maximum accuracy.

3.6 Sugar correction

For samples with **high sugar content**, sugar correction must be applied to obtain correct results. High sugar content in samples lowers the concentration of methanol in the headspace. As a result the Spark M-30 will overestimate the methanol concentrations of such samples.*

Use one of the following approaches to counter this effect and obtain accurate results:

1. **Set the sugar content** of the samples in the settings menu (Settings → Advanced → Sugar correction). Measurement results will be automatically corrected, and the correct methanol and sugar concentration displayed in the result screen.
2. **Import results** into the Spark M-30 Data Spreadsheet (available on www.alivion.ch) and manually input the sugar content of your measured samples.
3. **Apply formula** provided below to correct results displayed by the device yourselves. Note that this requires that measurements were performed with a sugar content of 0 g/L set in the settings menu.

$$M = M_0 \cdot \frac{E \cdot 0.0114 + 1}{E \cdot 0.0114 + \frac{S}{1000} \cdot 1.95 + 1}$$

M: Corrected methanol concentration, ppm
M₀: Methanol concentration shown on the detector, ppm
E: Ethanol concentration of the sample, vol%
S: Sugar content of the sample, g/L

* For instance, the detector will show 240 ppm instead of 200 ppm for a sample with 0 vol% ethanol and a sugar content of 100 g/L.

3.7 Warnings

The Alivion Spark M-30 continuously monitors the environment and corrects the results for temperature and ambient air variations to ensure accurate measurements.

If maximum accuracy cannot be guaranteed, the following warnings will be displayed.

Calibration warning

Cal

A new calibration should be performed. Displayed if the current calibration is older than **24 hours** or if the calibration temperature differs from the ambient temperature by more than **5 °C**, or when the last calibration was performed with an active temperature or contamination warning.

Temperature warning

Tem

The temperature is not constant. Make sure you are in a temperature stable environment and wait for the device to equilibrate.

Contamination warning

Con

High background concentration in the ambient air. Measurement accuracy is decreased at low methanol concentrations. Remove solvents from the environment, ventilate or switch location. After, the device can be cleaned from impurities by starting a cleaning procedure from the settings menu (see chapter 5 for details).

WARNING

Performing calibrations and measurements with active warning flags is discouraged since imprecise results must be expected.

4 Data

- Measurements are saved in **measurement series** of maximum **100 measurements**.
- Up to **10 measurement series** can be saved on the device and exported to a PC.
- Only measurements from the active (most recent) measurement list are displayed on the device.
- Measurements can only be saved to and deleted from the active measurement list.

NOTICE

If the active measurement series is full, or to organize your measurements (by date or type), a new series can be created in the settings menu (chapter 5).

4.1 On the device

Saved measurements can be accessed in the *Data* menu accessible from the home screen:

- Measurements are listed in chronological order.
- A detail page of each measurement can be accessed from the data list via the confirm button.
- Individual measurements can be deleted from the detail page.

4.2 Export

Measurement series are saved as CSV-Files on the device and can be exported to a PC:

1. Connect the Alivion Spark M-30 to the USB port of a computer using the provided USB cable.
2. Copy and paste the data files from the device folder (in your file explorer) to another location.

**NOTICE**

Each row in a data file contains a measurement with values separated by commas.

**NOTICE**

The data files cannot be edited or deleted from the PC. Delete measurements individually from the *Data* menu (see chapter 4.1) or collectively via the *Resets* menu (see chapter 5).

5 Settings

To modify the *Settings*, select the corresponding item in the main menu.

Language

The following languages are supported:

English, German, French, Spanish, Italian

Date & Time

Set date & time.

Brightness

Set screen brightness. Lower screen brightness will increase battery time.

Key sounds

Activate or deactivate the sound when pressing keys.

Advanced

<i>Clean sensor</i>	Cleans the sensor after a contamination warning (lasts 15 min).
<i>Sugar correction</i>	Set the sugar concentration (g/L) of your samples for automatic correction of measurement results (see Chapter 3.6).

Resets

<i>New measurement series</i>	Creates a new measurement series starting with ID 1.
<i>Delete measurements</i>	Deletes all saved measurements from the device.
<i>Reset calibration</i>	Delete the current calibration and apply the factory calibration.
<i>Restore factory settings</i>	Delete all measurements and reset the calibration.

Software update

Used to update the device firmware. Firmware updates and instructions are supplied from www.alivion.ch.

About

Displays device information.

6 Device faults

An error message will appear if a device fault occurs. Please see chapter 7 (Troubleshooting) instructions.

If the fault persists even after the device has been switched on and off several times and you followed the instructions in chapter 7, contact Alivion's technical support.

7 Troubleshooting

This section explains device faults and provides instructions on how to resolve them.

7.1 Sampling errors

Timeout during sampling

- Remove the sample vial within 10 seconds after the Alivion Spark M-30 displays *Remove sample*.
- Make sure to confirm the removal of the sample vial within 10 seconds.

Sampling vial not connected correctly

- Make sure an Alivion sample vial is screwed into the Alivion Spark M-30 before pressing *Continue* on the *Insert sample* screen.
- Make sure only original Alivion sample vials are used.
- Make sure the vial is screwed in hand-tight.
- Retry the measurement with a new sample vial.
- Check if the sealing O-ring (see chapter 2.1) is missing. If so, replace it with the included replacement or otherwise contact Alivion's technical support.
- Check if the windings are damaged in any way. If so, please contact Alivion's technical support.
- Contact the Alivion technical support if the problem persists

Sampling failed (ERR 942)

- Do not remove the sample vial before prompted to do so.
- Ensure that the sample vial is not moved or tilted during sampling.
- Make sure the vial is screwed firmly and try using a new vial.
- Make sure the sealing O-Ring is installed correctly.

- Contact Alivion's technical support if the problem persists.

Sampling failed (ERR 943)

- Make sure only Alivion vials are used and try using a new vial.
- Contact Alivion's technical support if the problem persists.

7.2 Measurement errors

Measurement failed (ERR 950/951)

- Ensure the inlet of the Alivion Spark M-30 is clean. Remove debris or liquid sample at the inlet of the device.
- Do not screw in a sample vial or block the inlet in any other way during the analysis.
- Contact Alivion's technical support if the problem persists.

Measurement failed (ERR 960/961)

- Retry the measurement with a freshly prepared sample.
- Retry the measurement with a different sample liquid.
- Try to clean the sensor (Settings → Advanced → Clean sensor) and re-calibrate the device.
- Contact Alivion's technical support if the problem persists

Critical error (952/953)

- Restart the device and retry the measurement.
- Contact Alivion's technical support if the problem persists.

Ethanol concentration too high (954)

- The Spark M-30 is only able to measure samples with an ethanol concentration up to 20 vol%.
- Contact Alivion's technical support if the problem persists.

7.3 Calibration errors

Calibration failed. Measurement error

- The calibration failed due to a measurement error described in chapter 7.2. Follow the instructions for the respecting measurement error.

Temperature difference during calibration too large

- Ensure the ambient temperature does not change during the first and second calibration measurement. Retry the calibration once the temperature is constant.

Calibration failed (ERR 972-976)

- Ensure the correct standard was used for the calibration.
For calibration without ethanol *Alivion Methanol Standard 1 WITHOUT ETHANOL* (orange liquid, red bottle cap) is used first, followed by *Alivion Methanol Standard 2 WITHOUT ETHANOL* (blue liquid, red bottle cap).
For a calibration with ethanol *Alivion Methanol Standard 1 WITH ETHANOL* (orange liquid, blue bottle cap) is used first, followed by *Alivion Methanol Standard 2 WITH ETHANOL* (blue liquid, blue bottle cap).
- Ensure only original Alivion calibration standards and sample vials are used.
- Ensure the inlet of the device is clean and free of debris or liquid sample.
- Ensure any used sample vials and pipettes are not near the inlet of the Spark M-30.
- Ensure no sample liquid is on the windings or rim of the sample vials.
- Try to clean the sensor (Settings → Advanced → Clean sensor), reset the calibration to factory settings (Settings → Resets → Reset calibration) and retry to calibrate the device.
- Contact the Alivion technical support if the problem persists.

7.4 General errors

Critical error (ERR 900)

- Restart the device and retry the measurement.
- Contact the Alivion technical support if the problem persists.

Charging not possible during measurements

- Remove the charging cable from the device and retry your measurement.
- Disconnect the Alivion Spark M-30 from the computer and retry your measurement.
- Otherwise restart your Alivion Spark M-30 without a cable connection.

Please clear inlet

- Make sure nothing blocks the inlet of the Alivion Spark M-30. Remove any sample vials or debris at the inlet.
- Contact the Alivion technical support if the problem persists

Temperature not within operational range

- The device can only be used for measurements within the operational temperature range of 10–35 °C.
- Make sure the ambient temperature is within the operational temperature range, let the device equilibrate and retry the measurement.

Measurement series is full

- A maximum of 100 measurements can be saved in the active data series.
- Create a new data series under:
Settings → Resets → New measurement series.

My results seem not accurate/consistent

- Ensure there are no warnings are active (see chapter 3.7).

- Ensure the device inlet is not contaminated with sample liquid. Otherwise clean it with a damp cloth, let it dry and try again.
- Ensure there is no sample liquid on the windings or the rim of the sample vials.
- Ensure the device, vials, pipettes and your sample liquid are all at ambient temperature.
- Ensure you prepare the sample vial with exactly 2 mL of sample liquid.
- Perform a new calibration.
- Try to clean the sensor (Settings → Advanced → Clean sensor).
- Contact the Alivion technical support if the problem persists

Contamination warning does not go away

- Make sure the device is operated in a room with clean air.
- Ventilate the room. If necessary, remove any potential contamination sources (e.g., open sample liquid, used sample vials and pipettes).
- Remove the O-ring from the device. Clean it and the sample vial inlet with a damp cloth. Dry both and put the sealing O-ring back in.
- Try to clean the sensor (Settings → Advanced → Clean sensor). Try using it repeatedly if necessary.
- Contact the Alivion technical support if the problem persists.

8 Maintenance

8.1 Transport and storage

Transport the device so that it is not exposed to shocks, moisture or direct sunlight.

Store the device in a dry, well-ventilated environment. Protect the device from direct sunlight and other heat sources.

In order to ensure the proper functioning of the device, it may only be stored under the following conditions:

Storage temperature:	5–30 °C
Relative humidity during storage:	max. 80 %

**WARNING**

Do not store the device in an environment containing alcohol or other solvents.

8.2 Cleaning

**WARNING**

Damage to the device! Never use solvents or solvent-containing cleaning products for cleaning. Do not immerse the device in liquids and do not allow liquids to reach the connections.

Turn off the device before cleaning and make sure it is not connected to any electrical circuits.

Clean the device with a dry, soft cloth. For heavier soiling, use only a very slightly damp cloth.

Moisture must not be allowed to enter the device. If moisture has nevertheless penetrated the device, do not use the device any longer.

8.3 Disposal

This product must not be disposed of as municipal waste. It is therefore marked with the adjacent symbol. Alivion takes this product back free of charge. You can get more information from Alivion and the national sales organizations.

Batteries and rechargeable batteries must not be disposed of as municipal waste. They are therefore marked with the adjacent symbol. Dispose of batteries and rechargeable batteries at battery collection points in accordance with the applicable regulations.

9 Technical data

Alivion Spark M-30

Measuring range	0-999 ppm* methanol
Accuracy**	5 ppm or 5% of the measured value (whichever is larger)
Resolution	1 ppm
Limit of detection	1 ppm
Ethanol Limit	0-20 vol% ethanol in sample
Operational temperature	10 °C to 35 °C (50 °F to 95 °F) Automatic temperature correction
Operational humidity	0-90 % relative humidity
Measuring principle	Headspace gas chromatography
Measurement time	2-4 min
Sample volume	2 mL
Internal storage	Up to 1000 measurements
Power supply	Rechargeable LiPol Battery (3900mAh)
Dimensions	195 x 100 x 40 mm
Weight	ca. 360 g
Interfaces	USB-A cable to PC
Protection class	IP31
CE - Certification	EMV, RoHs, UN dot 38.3

* The device informs the user when the concentration exceeds 999 ppm.
Mass concentration (DIN 1310) at a reference temperature of 25 °C.

** Standard deviation of the reproducibility with methanol-water mixtures (10-1000 ppm).

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The screenshot shows the measurement interface of the Alivion Spark M-30. At the top, the time is 14:35 and the battery level is indicated. The main display shows the following data:

Nr.	Datum	Zeit
24	19.03.2022	15:01

Below this, the primary measurement is displayed: **Methanol, ppm** with a large value of **240**. To the right, the temperature is **Temp, °C** at **25.2**. Below the temperature, the sugar concentration is **Sug, g/L** at **150**. A warning icon labeled **Tem** is shown. A **Kalibration** section is visible below the warning, with its own **Datum** (19.03.2022), **Zeit** (14:15), and **Temp, °C** (25.0). At the bottom, there are two buttons: **LÖSCHEN** (trash icon) and **SPEICHERN** (document icon).

Annotations on the right side of the screenshot point to the following elements:

- Details zur Messung (points to the top measurement table)
- Methanolkonzentration (points to the 240 ppm value)
- Temperatur (points to the 25.2 °C value)
- Zuckerkonzentration (Kapitel 3.6) (points to the 150 g/L value)
- Warnungen (Kapitel 3.7) (points to the Tem warning icon)
- Details der entsprechenden Kalibrierung (points to the calibration table)

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All rights reserved. Edition October 2023
Alivion AG, Technoparkstrasse 1, CH-8005 Zürich